

FUJIIRYŌKI

*Comfort to your satisfaction*



**CYBER-RELAX**  
Massage Chair SKS-1650

# Embraced by pleasantness on a sunny afternoon

If possible, we should not carry over the day's fatigue to the next day.

With the Solution massage, the whole body is wrapped up pleasantly.

It sets free the fatigue felt at the end of day.

While gently kneading down the whole body, firm stimulation is given to areas where stiffness is felt.

Enjoy a luxurious time to heal fatigue experienced in everyday life.

Toward liberation from stiffness of muscles and fatigue



heal your body, release yourself **CYBER-RELAX**



*The happiness  
of being embraced*

## Gentle relieving of stiffness from hands and arms

### **3D mode massaging**

Employment of 3D mode massaging that incorporates movement of the entire arms in addition to the hands. Variegated and 3D kneading courses that incorporate rear and front movements with such basic functions as wavelet, tapping, knead up and knead down reproduce a sensation previously offered only by expert masseurs.

### **Variegated kneading skills and 509 types of massage**

3-D mode that adds up rear and front movements of knead balls to the basic capacities of air, knead, tapping, stretch, shiatsu massage, etc. If the rolling is combined, it provides up to 509 types of massage.



### **Arm massaging**

Grasps, presses and kneads. Gentle stimulation of air bags produces a pleasant feeling while peacefully relieving the fatigue accumulated in the hands and arms.



By changing reclining angles, a wide area from arms to hands can be massaged.



*Flooding Joy*

Great care at the lower half of body where fatigue is likely to accumulate

**Seat surface massaging**

Soft stimulation over a surface as if being pushed by air relieves stiffness from the muscles of the hip to the back of calves, and promotes the blood circulation around the ankles.




**Sole massaging**

Stimulates the soles where many massage points are concentrated, and promotes blood circulation by relieving stiffness from the legs where fatigue of the day tends to accumulate.



**Calf massaging**

Stimulation resembling a milking action enhances the “milking” function that helps return the blood that tends to stagnate in the lower half of the body.



**New foot stretch**

Stretching of the knees and the lower half of body, while grasping and holding the legs, pushes the calves up and down slowly with the air bag on the seat surface. It massages the knee joints with movements to stretch them gently, to promote blood circulation in the legs.






Various kinds of information are displayed on the screen in easily comprehended ways.

**Air Massage Selector button**

**ARM MASSAGER**

**Pulse Mode**

Pulse Mode combines the gentle massaging with air with delicate stimulation just like by hands. Blood circulation is accelerated with rhythmic stimulation of 7 cycles comprising pressurizing, sustained pressurizing and release of pressure.



Open the cover

**Manual control panel**



\*All illustrations are images.

Fine interior article  
to suit your taste

# CYBER-RELAX

Massage Chair SKS-1650



### ■ Casters for moving the location

The chair can be relocated easily.



### ■ Removable back and seat pads, and pillow

They can be removed easily for maintenance.



### ■ Sole Shiatsu pad (Accessory)

Antibacterial specification sole pads for use on the bottom of the footrest.



AS-TS2

### ■ Side table (Optional)

A drink, newspaper or some small article can be kept within reach.

## ⚠ SAFETY PRECAUTIONS

■ Read the USER MANUAL carefully beforehand to use the product properly.

■ The following persons should not use the product, as it could cause physical problems. ● Persons who are prohibited to receive a massage, for example, those with thrombosis (embolism), severe aneurysm, acute varicosity, various dermatitis and skin infections (including inflammations of subcutaneous tissue), etc.

■ The following persons must consult a doctor before using the product, as its use could cause physical problems. ● Those who suffer from malignant tumors. ● Those who are in an unstable condition in early pregnancy or immediately after childbirth. ● Those who suffer from anesthesia related to the peripheral circulation trouble due to diabetes, etc. ● Those who have a body temperature higher than 38°C (pyrogenic period) (Example: When strong acute inflammatory symptoms are observed, such as malaise, chill, anisopiasis, etc. When the person feels weak.) ● Those with wounds in the locations to be treated. ● Those who suffer from acute illness (dolorific) ● Those who suffer from the osteoporosis. ● Those who suffer from inflammatory diseases such as sprains, muscle strain, etc. ● Those who suffer from a sensory disorder. ● Those who have a heart disease. ● Those who use a medical electronic device embedded in the body such as a pace-maker, which is susceptible to electromagnetic interference. ● Those who have some abnormality in the spine or whose back is bent to the right or left. ● Those who suffer from anemia. ● Those who need to massage a part of body treated medically in the past. ● Those who need repose. ● Those who feel any abnormality of the body other than the above.

### ■ Motor-driven reclining

Angles of the backrest and footrest for the legs can be adjusted with the remote controller from the sitting position.



### Specifications

Measurement (approx.)	(Footrest stored) Width 76 × Depth 123 × height 120 (cm)
	(Top of footrest reclined) Width 76 × Depth 196 × Height 78 (cm)
Weight (approx.)	70kg
Reclining angle (approx.)	120 - 170°
Power supply	(50/60 Hz)
Rated time	30 minutes
Power consumption	Main unit - 120 W, heater - 20 W (50/60 Hz)
Medical Equipment Approval in Japan	27BZ0878

Fuji Medical Instruments Mfg.Co.,Ltd.

● For further information;

\* Design and specifications are subject to change without prior notice for improvement of the product.

\* Colors of printing and actual product may differ slightly.