

Select programs with the easy-to-use control.



There are 4 massage programs for you to choose from.

- 1 **Full-body Massage & Stretch**
Full-body massage & stretching procedure is used.
- 2 **Stretch All**
A procedure for full-body stretching for the building of correct postures.
- 3 **Full-body Massage**
A procedure which can help listless and easily fatigued people regain their vigor.
- 4 **Lower Body Massage Original**
A procedure that focuses on the waist and lower body. Recommended use: get rid of muscle fatigue below the waist to around the pelvis as well as leg fatigue.

Warm-up
Used before exercise to flex muscles or in the morning to wake up the whole body.

Cool-down
Used after exercises to balance the body, relax muscles or stretch, so as to avoid post-workout fatigue the next day.

Max intensity for a more powerful massage with the touch of a button.

For Ultimate "Convenience of Use"

■ **Headrest Cover** Made of easy to clean materials.

■ **Extendable leg rest** Can be freely stretched or contracted according to body height.

■ **Installation at Home** Minimum width of the chair: 55cm *When the headrest is taken off.

It can be moved in when the passage is narrower than 65cm.

Height of Headrest About 55cm when the headrest is taken off.

Set the chair upright and move it in. About 55cm.

About 163cm

About 61cm

* Dimensions at the time of unpacking
* Install the armrest after the chair is moved in

Specifications



The length is about 135cm and the weight is about 74kg.

| | |
|----------------------|--|
| Product Name | i.12 |
| Power Supply | AC220-240V ~ |
| Rated Power | 150W |
| Outer Layer Material | Synthetic Leather |
| Angle of Inclination | 115°-165° to the level surface; |
| Accessories | AC power cord, acupressure sole pads (L/R), 2 acupressure plates, remote control bag, L-shaped spanner |

Colour



The History of Massage Chairs is also that of Inada.

Which Has Been a Professional Manufacturer of Massage Chairs for Half a Century.



Based on preventive medicine, we bring health to people all over the world through development of massage chairs

- Actual colours may differ slightly from colours in the catalogue due to printing.
- Please keep in mind that specifications are subject to change without prior notice due to product improvements.
- The included information may differ from the specifications of this product.
- Although the term "shiatsu" is used on the remote control and accessories for this massage chair, no alternative for shiatsu is included in the efficacy or effect.

Member of The Japan Home-Health Apparatus Industrial Association



MEDICALCHAIR
Stretch & Smooth for S-Line

i.12



Oakdale

Official European Importer/Distributor

Ellesmere Business Park
Oswestry Road, Ellesmere SY12 0EW, UK
Tel: +44(0) 1691 623234 Fax: +44(0) 1691 624224
Freephone: 0808 1682231

Web: www.inada.co.uk Email: sales@inada.co.uk


* Product information in this catalogue is as of September 2016

"Better Than Ever" Comfort and Quality of Massage



New Roller with "12 Acupoints"

The new "12 Acupoint" rollers are designed to strengthen the overall feeling of acupressure. The acupoints are the same size as the actual thumb.




Quality in detail


The chair is designed with comfort in mind. Inada pursues a philosophy of "quality in detail", therefore a mere touch is an excellent representation of its maximum comfort.



Equipped with "Seat Surface Airbag Functionality"



A balanced and therapeutic massage is achieved through a squeezing action, which is felt between the surface airbag and the heart.



When the chair is reclined a combination of massage rollers and reclining motion enables an entire body stretch.



Equipped with heaters for fingers and toes to warm extremities

The body is relaxed from the inside by gently warming fingers and toes that always feel like the coldest parts of the body.

Stretching - Continued improvement of Inada technology

Equipped with Full Stretch Function



Stretches the inner abdominal muscles, and relaxes the diaphragm

Equipped with Posture Correction Function



Firmly stretches the back by arching it entirely

Accurately detecting the acupoints for massage that have an effect on health

The smart optical detecting system can sense acupoints according to body shape

As the most important feature of the INADA medical massage chair, the smart optical detecting system locates acupoints and treats them.

Locating the acupoints in this way enables the chair to provide a specific massage according to the body shape



● Trapezius Massage

Deep massage of Trapezius


The chair massages along the trapezius, while avoiding the scapula, enabling treatment to parts of the trapezius which is prone to soreness and fatigue.



● Relaxing the Neck

Relaxing the muscles around the Tain zhu point through neck massage


The chair relaxes the Tain zhu and Feng chi (Wind Pool) points by using effective technique to alleviate stiffness and soreness.



● Neck and Shoulder Massage

Precise massage


The chair provides deep massage of the shoulder well where fatigue easily accumulates. Precise massage will alleviate tension and soreness around the massage well.



● Pelvis Massage

Squeezing the pelvis to relax muscles

The chair provides a wrapping massage from the waist down to under the hip, eradicating aches and tenderness around the pelvis through effective massage movement.



● Deep Massage of Lower Back

Deep massage of the lower back to relieve fatigue

The chair gives a long range of treatment, deeply massaging the muscles around the hip and relaxing every muscle between the waist and the pelvis. The lateral range of treatment is a wide as 20cm. This wide range provides an effective massage of the outer parts of the waist relieving fatigue.



● Massage of Legs and Feet

Increased blood circulation through air pressure from three directions

The chair's airbags which cover the legs from the calves to the soles of the feet, can ease soreness. At the same time acupressure pads are used to stimulate the soles of the feet so as to provide a strong and deep massage.



Acupressure Pads for Soles of the Feet



● Massage of Hands and Arms

A gentle but thorough chest spreading movement whilst relaxing the arms


From the elbow to the fingertips the chair performs a gentle but thorough air massage while comfortably stretching the chest soothing aches and tension in the arms.



Acupressure Pads for Arms

■ Wrapping Style Massage

Firmly gripping the body from shoulder to arm for a deep tissue massage



The chair's airbags are able to wrap around the body providing a deeper massage.

■ Professional Massage Theory


Automation of the massage programmed using professional massage theory



Professional masseurs will schedule treatment using massage theory according to each symptom and perform detailed massage procedures to each person according to their symptoms.

■ Four Roller See-Saw with Ultra Slow System

Recreating the feeling of double hand pressure by a professional masseur



Developed to give a comfortable and stimulating roller massage by breaking down body pressure to four points. Slow movements with brief pauses imitate those of a professional masseur and can relieve severe fatigue.