

# **FMI-EUROPE AB**

European representative of

*Fuji Medical Instruments MFG. CO., LTD.*

## **CYBER-RELAX**



*The Pioneer in massage chairs since 1954*

The standard model in a stylish design  
equipped with advanced functions

# EC-2700

## Unique Function

Making a diverse and sophisticated  
air-massage possible

## エアマジック Air Magic

The high-capacity, high-output  
compressor effortlessly controls  
35 air bags to provide Fujiiryoki's  
unique and innovative massage.  
The whole body is soothingly  
massaged by complex actions  
such as stretch and twist.



Combination of  
mechanism and air **677** types of massage

The combination of Fujiiryoki's Super-Mechanism with two massage balls and air massage makes it possible to provide 677 types of fine massages from the nape to the soles.

## Unique Function

## KIWAMIMECHA Massage



### Neck massage

The massage balls  
improve blood circulation  
around the neck by  
gripping and stretching  
the nape.




### Super Knead / Super Tapping

The tightness in the  
trapezius is relieved  
by providing deep  
stimulation to the upper  
part of the shoulders.



### Super Back Knead / Super Back Tapping / Buttocks Massage

These functions are  
effective in soothing  
the soreness caused  
by the strain in the back  
or buttocks, neuralgia  
and muscle ache.



### Firm grip and slow stretch

#### Stretch Massage

The calves are massaged by the air bags in the ottoman by gripping and slowly stretching them. This massage improves the circulation of the blood in the lower body, which is prone to swelling.



#### ● Stretch

The legs are held firmly and stretched forward.



#### ● Leg Stretch

The legs are held firmly and stretched downward.

#### ● Double Leg Stretch

The legs are held firmly, stretched forward and then downward.

### The lower back is sterically stretched by the cross & twist action

#### Twist Massage

The air bags located in the area of the chair for the lower back and buttocks independently move from side to side, up and down to provide a 3D massage. It loosens tense muscles in the lower back and buttocks.



#### ● Air Cross

Air to the lower back and buttocks is circulated diagonally to provide a firm stretching action.



#### ● Air Twist

Air to the lower back and buttocks on one side of the body is circulated simultaneously to provide a gentle stretching action.

### Hand & Arm Massage

There is a concentration of pressure points in the hands and arms. Tension in these points is properly released by two types of massage: Grip and Knead.



### Leg Massage

The air bags aid the "milking action", which promotes the circulation of blood in the lower body, by repeating compression and release from the lower part of the calves to the upper part of the calves, much like milking a cow.



### Buttocks Air Massage

Soft air stimulation gently loosens the muscles in the buttocks and the back of the thighs to promote a better blood circulation in the groins.



### Sole Massage

The soles are stimulated by air bags imbedded with pressure balls as the feet are firmly held from the heels.



### Lower Back & Buttocks Massage

The interlocking action of intensive massage to the lower back and air massage to the buttocks makes it possible to provide relief to the lowermost part of the sidebone, which cannot be reached by ordinary massage.





**Remote Control**

Sophisticated functions can be operated simply by using the remote control with easy-to-understand user-friendly operating buttons.

**Easy Quick Start**

Quick Start a session by choosing your height from 3 different settings, or 'My Memory' in which you can register your course.

**Register your course in 'My Memory'**

Register a course individualized to you according to the shape of your back and the positions of your shoulders and lower back.

**Cross-shaped Key for a simple operation**

Choices can be easily made from a variety of functions by the use of a user-friendly cross-shaped key.

# CYBER-RELAX



EC-2700 BK (black)

**Electric Reclining**

Find a comfortable, desired position via the easy to use remote control.



**Accessories**

**Removable Back Pad & Pillow**

The pad and pillow can be removed for easy cleaning.



**Sole Pads**

The antibacterial sole pads can be attached to the base of the ottoman.



**Specifications**

●Measurements (Approx) : [When ottoman is stored] W750×D1,300×H1,170(mm) [In reclining position] W750×D2,040×H720(mm) ●Reclining angle (Approx) :120 ~170 degrees ●Weight (Approx) :81kg ●Voltage : AC110V,120V,220V,230V,240V(50/60Hz) ●Power consumption : 115W ●Rated time : 30 minutes EC-2700's Japanese model AS-730 is 220AGBZX00123000 certified for medical device (Managed-Care Device) in Japan.

Model Number / Exterior	Upholstery	Remote Control	3D Point Navigation System	3D Mode Massage	Extreme Massage	Extreme Hip Massage Alignment	Fluctuation Massage	Waist/Hip Massage	Twist Massage	Sliding Footrest	Perfect Fit Footrest	Double Stretch	Lower Body Stretch		Sole Pressure Balls	Air Massage				Auto-course	Massage Types	Mechra Massage Speed Adjustment	Width Adjustment	Air Massage Adjustment	Memory Function	Skip Function	Repeat Function	Pulse Mode	Rhythm Mode	Simultaneous Air Massage	Remote Control Stand	Arranges	Reclining	Dimensions (mm) (Reclined Position)	Weight	Power Consumption
													Foot Stretch	Stretch		Shoulder	Arm	Back / Waist	Hip / Thigh																	
EC-3700	Leather (Black / Beige)	 ●User-friendly display ●Guidance Function ●Cross Key	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	W750×D1240×H1210 <W750×D2010×H780> 120-170° (Approximate) -Inclining Chair back and footrest -Independently Adjustable footrest (up/down)	86 kg	115 W
EC-2700	Leather (Black)	 ●User-friendly display ●Guidance Function ●Cross Key	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	W750×D1300×H1170 <W750×D2040×H720> 120-170° (Approximate) -Inclining Chair back and footrest -Independently Adjustable footrest (up/down)	81 kg	115 W
EC-1700	Leather (Black)	 ●Guidance Function ●Cross Key	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	W710×D1180×H1200 <W710×D1780×H720> 120-170° (Approximate) -Inclining Chair back and footrest -Independently Adjustable footrest (up/down)	79 kg	115 W

# NEW MASSAGE CONCEPT

**MASSAGE**  
(Kiwami mecha, Air)

- \* Deep massage with KIWAMI MECHA massage heads.

**NEW CONCEPT**  
"Air Magic"

- \* Total Body Stretch
- \* Twist (Twist/Cross)
- \* Double Stretch
- \* New Shoulder massage
- \* New Arm massage

## MASSAGE + AIR MAGIC

# FMI-EUROPE AB

European representative of

*Fuji Medical Instruments MFG. CO., LTD.*

Strandvägen 40A, 444 41 Stenungsund, Sverige, Tel: +46-(0)31-61 080 63

Mail: [info@fmieu.com](mailto:info@fmieu.com)

Web: [www.fmi-europe.com](http://www.fmi-europe.com)

Retailer:

### HISTORY



1954



1965



1975



1987



1995



2000



2001



2002



2004



2008

*Fuji Medical Instruments*

# CYBER-RELAX

